VISCOSITY SENSORS 440 & 443





Viscosity Sensors 440 & 443

The Viscosity Sensors 440 and 443 installs into a high pressure pump through type system. Both ends mate to HiP HF4 fittings with a 9/16-28 thread. The Viscosity Sensor 440 is rated for continuous operation up to 20,000 psi while the 443 sensor is rated up to 30,000 psi.

Fluid flows through ¼ inch high pressure port for measurement. Fluid is typically measured under static

conditions. A continuous flow adapter is optionally available. System measures viscosity, temperature and temperature compensated viscosity. It provides data via RS-232.

The high-pressure viscosity sensors are used in various applications including PVT oil research, EOR, core analysis, and supercritical fluid research.

Our industry leading viscosity technology

is based on a simple and reliable electromagnetic concept. Two coils move the piston back and forth magnetically at a constant force. Proprietary circuitry analyzes the piston's two-way travel time to measure absolute viscosity. A built-in temperature detector senses the temperature of the sampling chamber.

KEY FEATURES

Proprietary and Proven Technology:

- Accurate and reliable requiring minimal sample size (less than 2 ml)
- · Rugged stainless steel design that is highly durable
- Mercury free
- 440 sensor is pressure rated to 20,000 psi; the 443 sensor up to 30,000 psi
- · Measures gas and gas condensates in addition to fluid sample

Easy to Use:

- Self-cleaning
- · One moving part, the piston

APPLICATION RANGE

Viscosity measurements from 0.2 To 20,000 centipoise and 0 To 20,000 psi; extended option available for up to 30,000 psi

STANDARD METHODS

Correlates with ASTM D445
Correlates with ASTM D7483




solidpartners provensolutions

U.S.A.

PAC, LP | 101 Station Landing| Medford, MA 02155 T: +1 800.554.4639 | O: +1 781.393.6500 | F: +1 781.393.6515 sales.usa-boston@paclp.com

PAC Authorized Representatives are also located in most countries worldwide. For more information visit www.paclp.com.



Overall Viscosity	0.02 to 10,000 centipoise (cP)
Viscosity Ranges	0.022cP, 0.2-2cP, 0.25-5cP,
	.5-10cP, 1-20cP, 2.5-50cP,
	5-100cP, 10-200cP, 25-500cP,
	50-1,000cP, 100-2,000cP,
	250-5,000cP, 500-10,000cP
Viscosity Accuracy*	±1.0% of full scale
Repeatability	±0.8% of Reading
Temperature Sensor	PT100
Wetted Materials	Inconel 718 / 17-4 Stainless Steel
Maximum Particle Size	25-360 Microns
Max Ambient Temperature	375°F (190°C)
•	
Max Operating Pressure	SPL440 - 20,000 psi (1,406 bar) / 19,000 psi (CRN)
	SPL443 - High-pressure option - 30,000 psi
Power	Requires VISCOpro 2000 electronics
Certifications	CE, Alberta Boilers Safety Association, CRN #0F6429.2

^{* ±} of 5% of full scale 500-10,000 and 1.000 to 20,000cP Ranges

OPTIONS

Extra High Pressure option- rated up to 43,000 PSI (3000bar)

Extra High Temperature - for high pressure viscosity analysis at temperatures up to 300° C

Multiple ranges, easy change over

Continuous flow adapter

Hastelloy or silicon coated pistons

Integrates directly into PumpWorks®, Labview®, and other laboratory software information management software products





Copyright 2013/2 PAC L.P. All rights reserved 00.00.209

Cambridge Viscosity

Cambridge Viscosity With more than 10,000 installations worldwide,

Cambridge Viscosity is the proven leader in viscosity management technology. With over 25 years of experience, Cambridge Viscosity understands and meets the needs of laboratory researchers and process engineers in a wide range of industries whose jobs depend on the quality, accuracy, and reliability of viscosity measurement equipment. With their patented sensor technology, Cambridge Viscosity has become the gold standard in small sample viscosity measurement.









