VISCOlab 3000





VISCOlab 3000

The VISCOlab 3000 contains a pistonstyle electromagnetic sensor and integrated thermometer that provides continuous viscosity and temperature readings.

Measurements can be made in any of twelve 20:1 viscosity ranges. Viscosity can be displayed in either centipoise or centistokes, temperature can be displayed in °C or °F. Integrated temperature control allows user defined measurements at any setting from above ambient to 356°F

(180°C) and from above ambient to -40°C with optional cooling jacket. The sensor is easily field calibrated. Multiple operating ranges and product set-ups can be precalibrated.

Data can be output to a PC, via the RS-232 serial port for tracking, storage, or graphical display. The electronics allow the user to define temperature and viscosity display units, data averaging interval, density and much more. The electronics also provide digital readings,

as well as statistics on sample stability that allow greater repeatibility between tests and operators.

The VISCOlab 3000 is available in application specific configurations for diesel, lube oil, bunker oil, asphalt/bitumen, and gear oil. See specifications chart for more information.

KEY FEATURES

Easy to Use and Maintain

- Digital displays of centipoise and centistokes
- Stores up to 40 product specific test set-ups
- Removable liquid recirculating jacket and cap

Enhanced Features

- Minimizes product waste by requiring less than 2ml of fluid
- Requires small amount of bench space
- · Lightweight and portable
- Internal heater for temperature control at 40/50/100/135/150

High Accuracy and Reliability

 Software measures the sample multiple times until the standard deviation is less than 0.2% before displaying a viscosity value per ASTM 7483

APPLICATION RANGE

Viscosity measurements from 0.2 To 20,000 centipoise

STANDARD METHODS

Complies with ASTM D7483
Correlates with ASTM D445



solidpartners provensolutions

PAC, LP | 101 Station Landing| Medford, MA 02155 T: +1 800.554.4639 | O: +1 781.393.6500 | F: +1 781.393.6515 sales.usa-boston@paclp.com

PAC Authorized Representatives are also located in most countries worldwide. For more information visit www.paclp.com.

SPECIFICATIONS

Viscosity	0.2 to 20,000 centipoise (cP)	
Viscosity Accuracy*	±1.0% of full scale w/ViscoPro 2000	
Repeatability	±0.8% of Reading	
Temperature Sensor	PT100	
Wetted Materials	316L and 430 Stainless Steel	
Maximum Particle Size	25-360 Microns (Range Dependent)	
Temperature Control	Above ambient to 356°F (180°C)	
Standard Cable Length	3 Feet (0.9 meters)	
Power	100-240 VAC, 50-60Hz	
VISCOlab 3000 Diesel : Measures diesel oil from 0.5-10cP at 40°C		
VISCOlab 3000 Lube Oil: Measures lube oil from 1-20cP at 100°C and 25-500cP at 40°C		
VISCOlab 3000 Bunker: Measures bunker oil from 50-1000cP at 50°C		
VISCOlab 3000 Asphalt: Measures asphalt from 250-5000cP at 135°C and 150°C		
VISCOlab 3000 Gear Oil: Measures gear oil from 50-1000cP at 40°C		

^{* ±} of 5% of full scale 500-10,000 and 1.000 to 20,000cP Ranges

OPTIONS

Insulator Sleeve (For operating temperature above 80°C)		
Insulator Cap (For operating temperature above 80°C)		
Cooling Jacket (For re-circulation fluids)		
Low Temperature Cap (Up to 50°C)		
Optional Piston Ranges	0.2-2cP, 0.25-5cP, 0.5-10cP,	
	1-20cP, 2.5-50cP, 5-100cP,	
	10-200cP, 25-500cP, 50-1,000cP,	
	100-2,000cP, 250-5,000cP,	
	1500-10,000cP*, 1,000-20,000cP*	



Cambridge Viscosity

With more than 10,000 installations worldwide, Cambridge Viscosity is the proven leader in viscosity maanagement technology. With over 25 years of experience, Cambridge Viscosity understands and meets the needs of laboratory researchers and process engineers in a wide range of industries whose jobs depend on the quality, accuracy, and reliability of viscosity measurement equipment. With their patented sensor technology, Cambridge Viscosity has become the gold standard in small sample viscosity measurement.



Copyright 2013/1 PAC L.P. All rights reserved 00.00.210















