



Accurate Viscosity Measurement under Extreme Pressure & Temperature

*Gulf Coast Conference
October 17, 2012
William A. Small*

Agenda



- *Cambridge Viscosity & PAC History*
- *Why High Pressure/Temperature Viscosity measurements required*
- *Viscosity Importance for Reservoir Fluid Analysis*
- *Reservoir Fluid Analysis Overview & Challenges*
- *Cambridge Viscosity Oscillating Piston Viscometer Technology*
- *Operating Characteristics of VISCOlab PVT Viscometer*
- *Pressure & Temperature Standards*
- *Summary*

Cambridge Viscosity & PAC History

- *Founded 1984*
- *Located near Boston/Cambridge, MA*
- *The Leader in Small Sample Viscosity*
- *Developed 1st High Pressure/Temperature Viscometer for Reservoir fluids in 1990*
- *Proprietary Technology*
- *Blue Chip Customer Base*
- *Over 10,000 Systems Installed*
- *Worldwide Reach*
- *Compliance with ASTM D7483 & D445*
- *CVI partnered with PAC in 2012*



About PAC

- *Part of ROPER INDUSTRIES - (NYSE: ROP)*
- *Leading global provider of advanced analytical instruments for industrial laboratories and online process analytics*
- *> 400 employees and over 30,000 instruments installed*
- *11 PAC offices and > 50 distributors Worldwide*
- *> 50 Engineers in Six Countries*
- *Core technologies including:*
 - *Gas chromatography*
 - *Physical properties*
 - *Fuel performance*
 - *Elemental analysis*
 - *Process analytics*
- *Consists of leading, long-established product brands*



Where is HP/HT Viscosity

Measurement Required

Upstream

- *Exploration - Fluid Reservoir Analysis*
- *Production Operations*
- *Enhanced Oil Recovery using Polymers (EOR)*
- *Corrosion Inhibitors Analysis*
- *Brine Solutions & CO2 Analysis*

Other Applications


- *Lubricants at “contact bearing pressure” Analysis*
- *R&D applications*
- *Ionic liquids Analysis*

Viscosity Importance for Reservoir Fluid Analysis

Importance of viscosity is two fold:

- Determining well productivity using Darcy's law
(describes flow through a porous medium)
- Analyzing displacement in recovery
(i.e. – amount of oil recovered vs. water injected)

$$Q = \frac{-kA}{u} \frac{(P_b - P_a)}{L}$$



Reference: Kurt Schmidt of Schlumberger, Standards Workshop, 2012


Oil Reservoir Characteristics

Drive Economics



Geological characteristics of the rocks and fluids - **control** the performance of hydrocarbon reservoirs

- Involves numerous properties of the particular rocks and fluids
- Including interactions between the rocks and fluids
- High quality data are critical for:
 - Accurate predictions of reservoir performance for initial investments and analysis over time
 - Selection of optimal depletion strategy over time



The Four Most Critical Production Factors

- **Original Hydrocarbon in Place:** Volume of fluid in the rock formation
- **Production/Injection Rates:** Flow rates through porous rocks for the hydrocarbon fluids in the particular rock formation
- **Fractional Flow:** Rate that water flows through the reservoir rock/fluid combination
- **Recovery Factor:** Portion of the hydrocarbon that can be extracted given the rock and fluid characteristics



Sources:

- 1) Honarpour, MM, et al, *Rock/Fluid Characterization and their Integration – Implications on Reservoir Management*, 2006, SPE103358
- 2) Dindoruk, B., *“Reservoir Fluid Challenges: Viscosity”*, 2010, Cambridge Viscosity-Schlumberger Roundtable Workshop on High Pressure and High Temperature Viscosity Standards

Accurate Viscosity Data: is Vital

Key Reservoir Characteristics	Significant Factors				
	Fluid Viscosity	Rock Permeability	Oil Volume (Formation Volume Factor)	Oil and Water Saturation	Rock Porosity
Hydrocarbon Volume			<i>b</i>	S	<i>f</i>
Production - Injection Rates	<i>m</i>	<i>k</i>	<i>b</i>		
Fractional Flow of Water	<i>m</i>	<i>k</i>	<i>b</i>		
Recovery Factor - Potential	<i>m</i>	<i>k</i>		S	


Viscosity Uncertainties Have a Significant Impact:
 Typical Variances without viscosity data -
 > 10 % for Light Oil and
 > 15% for Heavy Oils

Source: Honarpour, MM, et al, SPE 103358

Reservoir Fluid Analysis Challenges



- Only **small volumes** of oil are extracted for evaluation
 - Typical “bomb” is 250 ml
- Material is **very expensive** to capture
 - Typical sample costs \$250,000
- Samples Often include:
 - **Corrosive** fluids (H₂S, brine, drilling muds & contaminants)
 - **Excessive Water** - leaving only 50 ml of oil to analyze
 - And are under - **High pressures & temperatures**



And For Analysis


- Fluid **sample conditions must be maintained** from reservoir to lab

Otherwise: gas flashes off, **irreversibly** changing fluid viscosity

Cambridge Viscosity Technology



“Electro-Magnetic Viscometer”



- *Piston motion completely controlled by two magnetics around sensor body*
- *Total Piston travel for 1 full cycle is measured in micro-seconds*
- *Calculation is made on the resistance of the piston in the fluid & time of travel of piston*
- *Viscosity is then calculated and presented in centipoise - cP*


“Sensor shown in a process line”

CVI VISCOLab PVT:

Ideal for Reservoir Fluid Analysis



- *Integrated PVT system:
Pressure – Viscosity – Temperature*
- *High pressures: up to 40,000 psi*
- *Small Sample size – only 6 ml*
- *High temperatures: up to 315 C*
- *Off-gassing incorporated*
- *Accurate: $\pm 1\%$*
- *Easy to use*
- *Manages corrosive fluids: H_2S , brine*
- *The only System that can measure Gas*



VISCOLab PVT has become the industry Gold Standard

VISCOlab PVT:




Single Integrated Interface

Traceable Pressure & Temperature
(PCV) Compensation

On Screen Statistical Analysis
of Test Conditions

Single Click Temperature Set

User configurable Dashboard &
Automatic Storage of Test Data




Easy to use.....

Pressure & Temperature Compensation

- **PVT Sensor calibration compensation is necessary**
 - Critical Sensor Feature: Annular space between piston & chamber wall
 - Space to increase with temperature and/or pressure
- **Compensation factor based on:**
 - Published data on fluids characterized by NIST and others
 - Data available to 20,000 psi and 150 degC
 - depending on fluid and test, but limited by viscosity, temperature and pressure

- **Compensation is based upon best-fit correlations for:**
 - Temperature and pressure are linear
 - Each individual piston range is compensated



Fluid standards for extended pressures & temperatures are needed

Sensor Enhancements & Options




- Easy piston change with **low torque** sealing
 - to extend the life of sensor
 - to prevent seal damage


- Standard measurement is in static conditions, **New low flow version** for continuous flow measurements is now available

- New H_2S –friendly materials (**Inconel, Hastelloy**)

- Extended Pressure & Temperature Capabilities are \longrightarrow coming....




Low torque sealing



Low flow version

Ultra-High Temperature and Pressure Systems



Now available


Key Points to remember On CVI PVT Viscometer

- **Widest** operational range with minimal sample volume used
- the only system that can measure gas
- Unique **automated** analysis
- **Simple** cleaning steps - increases lab efficiency
- **Robust**, Ergonomic Design Menus - lead to excellent analysis performance
- **Delivered** completely tested, calibrated and factory assembled for immediate turnkey operation
- **Very Cost Effective** – accurate results the first time



High Return on Investment

- Based on the **sample size** & cost per year - compared to other systems



- Based on amount of oil sample available after water is removed - **typically only 50 ml**
- Achieving accurate results with the first analysis

Users in the Oil Market today

- **Global** & Independent worldwide Oil Companies
- Oil Field Service Companies
- **Chemical** Manufacturers
- Fluid Reservoir Labs
- Oil Exploration Companies using Enhanced Oil Recovery (**EOR**)
- **Universities** & Scientific Institutes – performing analysis of phase change under extreme conditions for a wide variety of fluids & gasses




Industry Feedback and Trends

- Operation under **lower flow** rates
- Improved resistance to **H₂S & Brine** fluids
- Increased Pressure and Temperature required
- Improved ease of use and **data collection**
- Extend the Viscosity measurement of **Gas**
- Assist in establishing HP/HT **standards**



Enhanced Oil Recovery - EOR

- Oil Exploration Companies now are looking to modify the **OIL Temperature** for easier extraction
 - They are now using viscosity analysis to see what changes occur by manipulating the temperatures
- CVI Viscometers are now being used to **analyze** oil temperature manipulation




Summary:

Oil Exploration - Reservoir Analysis

VISCOlab PVT

– Industry Gold Standard:

- Accuracy
- Small sample size requirement
- High Pressure & temperature
- Ease of use
- Small footprint
- Wide viscosity range
- Fast analysis time
- Rugged System design
- Turnkey System
- Data Management



- Enhancements Continue

State of the Art -

High Pressure/Temperature Viscometer

Contact PAC

- For any questions, you can contact me at william.small@pacip.com. You may also email PACsales@pacip.com.
- Additional resources can be found at our website, www.pacip.com.
- Questions.....

Thank you



solidpartners
provensolutions