

## PAC Antek MultiTek Training Course

Office: +1.800.444.1803

Fax: +1.281.940.1803

pacsales@pacpl.com

### Course Description

This hands-on laboratory and lecture course covers a variety of topics to increase your understanding of PAC's MultiTek Elemental Analyzer. You will receive presentations for theory and component identification, mixed with practical hands-on training using the MultiTek instrument in our laboratory. Also covered will be the software for user interface, calibration, maintenance, and troubleshooting for common problems. Customized applications or specific peripheral combinations can also be accommodated upon request. These courses are conducted by certified skilled PAC technicians who have real world MultiTek Elemental Analyzer applications experience, much of which has been gained by visiting our customer base. You will also have access to the PAC applications and engineering teams for specific discussions and forums regarding your unique situations.

### Course Part Number and Cost:

Part Number: TR-MULTITEK

Cost: \$ 3,400 USD for each attendee (**INCLUDED: One Galaxy Tab A 10.1 with S Pen 16GB (Wi-Fi)**)

### Who Should Attend?

This class is intended for customers who have recently purchased a new MultiTek, new users, or current operators who need a refresher to fill in knowledge gaps.

### Length:

This class is 3 days

### Course Materials:

PAC will provide the presentations, as well as the instruments.

### Course Agenda:

#### Day 1

- MultiTek Configurations
- MultiTek NS Overview
- MultiTek NS Software Installation and Configuration

#### Day 2

- MultiTek NS Calibration and Sample Analysis
- MultiTek NS Data Analysis
- Discussion with Applications / Engineering

#### Day 3

- MultiTek Optimization
- Sample Handling
- MultiTek Troubleshooting

### Contact us to Register Today!

Web: visit our website [Online Registration Form](#)

Phone: 1-800-444-TEST or 281.940.1803

Email: [PACHOU.Service@pacpl.com](mailto:PACHOU.Service@pacpl.com)

### Terms and Conditions:

Payments need to be settled 100% prior to the start date of the training course.

Cancellation within 10 business days before the course starts or in case of non-participation, 100% of the course fees will be charged. Training registration can be transferred to a colleague.

Classes are limited to 8 participants. Courses with less than 2 registrations are subject to cancellation 10 business days before scheduled date.

Training price includes training materials and lunch/beverages at the training facility. Cost for travel, accommodations and other meals or expenses are not included.